

MIND IF WE SPEAK IT?

GENERAL i3 ADVANTAGE™ BENEFITS AND REQUIREMENTS

- ✓ We support interfaces with almost all EMRs, including Cerner, Meditech, McKesson, All Scripts, and Vista
- ✓ Solutions are scalable from single towers to large multi-facility enterprise systems
- ✓ Physical and/or Virtual Servers Supported
- ✓ High availability SQL, Server, and Storage solutions supported
- ✓ Multi-facility logic supports any number of sites on a single database
- ✓ Microsoft Best Practice. Security supports integration with Microsoft Active Directory
- ✓ Flexible user authentication options support for kiosk, 2 factor, CAC, and SSO
- ✓ Runs your anti-virus, connects to your domain, allows you to push Microsoft patches without calling us
- ✓ Physician report writer supports Windows Terminal Server, Citrix, or other virtual technology
- ✓ Windows 7 and 8 supported (x86 or x64)
- ✓ i3 Advantage has a Certificate of Networthiness (CON) for Department of Defense implementations
- ✓ Dell Data Protection Enterprise (DDPE) Encryption provided with our supplied workstations helps you stay compliant and protect PHI
- ✓ Dell 3-year on-site Pro Support provided with our supplied servers and workstations
- ✓ Microsoft Windows Server / Microsoft SQL Server Backend
- ✓ Xerox 3-year on-site support provided with our Xerox ColorQube printers
- ✓ ADT, SIU, and ORM for Inbound HL7 Demographics
- ✓ Flexible and customizable HL7 and other interface solutions supported
- ✓ ORU, MDM, and document cold feeds supported for results
- ✓ Remote VPN access required for system setup and ongoing support
- ✓ Windows fax server connected to dedicated analog line required for use of faxing interface option

INTERESTED IN IMPROVED COMPLIANCE AND A STRONGER BOTTOM LINE?

Find out if **Summit Imaging's i3 Advantage** is right for your Interventional Pain Management Team. Call **866-578-5713** or go to **summitimaging.com** to request a demo.



Summit Imaging, Inc.

306 SE 291 Highway, Suite 2
Lee's Summit, MO 64063

i 866.578.5713

e i3advantage@summitimaging.com

w summitimaging.com